

A REVOLUTION IN WELLNESS



THE TOTAL BIOLOGICAL TERRAIN RESEARCH STATION.

SEE THE PICTURE, GET THE NUMBERS.

\$9,995 pkg price
(Monitor sold separately)

Zooming Video Microscope

Unparalleled performance with superior imaging gives you the qualitative “look” at the underlying biological terrain.

Biomedx integrated fiberoptic microscope base with 150 watt quartz halogen lightsource, universal style turret condenser for phase contrast, darkfield, brightfield and 3D imaging, quintuple objective nosepiece with basic 4x & 40x non-oil phase contrast objectives, optical zoom system, binocular eyepiece assembly with trinocular port, deluxe rubber eyecup 10x FN 20 eyepiece oculars, multi gamma high resolution color sync DSP medical micro-head video camera.

BEV Multi-Parameter Analyzer +

A fast and simple process to see the “body-electric”. Gives you the quantitative “numbers” behind the picture - pH, redox, resistivity and a whole lot more – so much more we call it Flow Systems Auditing.

The Biomedx BEV bio-electronic testing apparatus comes with peripheral items - like a refractometer, urotensiometer, urinometer, reagent kits, training guide and real time telephone support. Merges the independent methodologies of Vincent, Reams, Revici and others into a powerful research tool so a clinician can quickly filter down to understand primary metabolic control points and focus on what matters most in the clients health.

BIOMEDX